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ENGLISH VERSION →


ATOMIC ABSORPTION SPECTROSCOPIE AAS

Atomic absorption spectrometry is a method of proportioning chemical elements founded on the absorption of atomic radiations in vapour phase.

MAIN CHARACTERISTICS

- Spectral range: 189 - 900nm
- 4 multi-elements lamps
- Band-width: 0.12nm with 2nm
- EDL function (*Electrodeless discharge lamp*)
- Automatic running of samples
- Automatic recognition of lamp type
- Deuterium reference lamp


OPERATIONS


The solution containing the element to be proportioned is aspirated through a capillary and directed towards the atomizer. It is then pulverized at very high speed through a room where the large drops are eliminated. The thus vaporized solution is pulled by fuel through the burner at the exit of which combustion takes place.

After thermal dissociation in the flame, incidental radiations provided by a hollow cathode lamp which emits the spectrum of the element that one wants to proportion are absorbed. The measurement of this absorption makes it possible to determine the concentration of the element to proportion.

REALIZATION EXAMPLES

Atomic flame absorption
Temperatures 3000°C

Spectrum absorption
Curve of absorbance according to the concentration

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